



The Noodle TA User Instructions

1 – **Never attempt to use The Noodle TA AFO without a shoe!!!** Put on a pair of socks. Loosen the shoestring and remove any shoe liners.



2 – Place The Noodle TA AFO into the shoe first.



3 – Replace the shoe liner on top of The Noodle AFO footplate.



4 – Place your foot into the shoe. Using a shoe horn may help.



5 – Slip the strap through the metal ring and fasten to the Velcro Strap-eezz pad.



6 – You have now properly donned The Noodle TA AFO.



Caring for The Noodle TA AFO

Washing – You may wash the AFO with soap and water. The cuff pad may be hand washed with soapy water and air dried. Do not use solvents, chlorine bleach or washing machines.

WARNINGS!!!

Never use The Noodle TA AFO without a proper shoe. Wearing the AFO without proper shoes can cause slips and falls which may result in injury.

The Noodle TA AFO should only be dispensed under the guidance and supervision of a qualified healthcare professional.